plasma ONE

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Table of Contents

1	Foreword	4
1.1	Intended use	4
2	Before using the device	5
2.1	Safety information	5
2.2	Checking the contents of the package	5
3	Commissioning	6-8
3.1	Commissioning and explanation of the energy control unit	
3.2	Commissioning and explanation of the converter	
3.3	Commissioning and explanation of the instrument probes	
3.4	Commissioning and explanation of the patient probe	
4	Warnings	9-12
4.1	Counter-indications	10
4.2	Cleaning, disinfection and sterilisation	11
4.3	Service & maintenance	11
4.4	List of plasma compatibility	12
5	Information for users	13-16
5.1	Successful treatment – step-by-step	
5.2	Experience in practice	14
6	Problem solutions	17
6.1	Solution table	17
7	Technical data	18-21
7.1	Explanation of symbols	18
7.2	Operating mode	18
7.3	EMC information	19
8	Certificates	22
8.1	ISO 13485 and ISO 9001	22
9	Warranty	23

Dear colleagues,
 Congratulations on your decision to include **plasma ONE** in your treatment concept.
 In part 1 you will find the plasma ONE instruction manual and in part 2 you will find practical information for users.
 We will be glad to answer any further questions about the product or its use by telephone, by e-mail or personally, if required. Your authorised dealer is also available to answer your queries if required. You can find the contract addresses of authorised dealers or the contact details of your personal medical product advisor on our Internet page.
 Please also consult the Internet page www.plasmamedicalsystems.com for further information about the latest studies and treatment methods, as well as further training.

This version does not claim to be complete and may be subject to technical amendments.

© 07/2013 plasma MEDICAL SYSTEMS® GmbH.

Reproduction of this documents or its contents is not permitted and required the explicit approval of plasma MEDICAL SYSTEMS® GmbH.

Version: 07/2013

1 Intended use

Dr. Susanne Pavisic

Dear colleagues, dear surgery team, to give you an overview of the future-oriented key technology which you have purchased, I would like to explain the intended use of plasma ONE.

The main idea behind our commercial philosophy and actions is the health of each individual member of the population as a whole.

The intended use of the medical product which you have purchased, plasma ONE, is disinfection!

Fields of application: dermatology, gynaecology, urology, surgery and dentistry.

User group: Physicians and specialist medical personnel.

The fact that plasma has antimicrobial properties is well known and has been scientifically researched and published throughout the world. The publications talk about **cold plasma** as well as **cold atmospheric air plasma**.

The causes of the antibacterial effect of plasma are due to heat, dehydration, shear forces, UV radiation, free radicals and charges. With cold plasmas, heat plays a subordinate role, as they will be practised at room temperature.

Particularly reactive radicals are generated in atmospheric plasmas. These include atomic oxygen, superoxide radicals, hydroxyl radicals, carbon monoxide and nitrogen dioxide.

These have a destructive effect on various cell components. If the cell walls of bacteria are directly exposed to the plasma, they become negatively charged due to bombardment by the electrons which are present in the plasma. The resulting electrostatic repulsion causes mechanical stresses, which can exceed the tensile strength of the cell wall and hence destroy it.

Bombardment with such particles can change or destroy the structural integrity of the cells. Cold plasma is the name which is given to plasma with a low temperature.

From a physical point of view, hot plasmas are in a state of thermal equilibrium with themselves and their environment, while cold plasma is not in such a state of equilibrium.

The plasma MEDICAL SYSTEMS® technology uses ambient air for ionisation.

A further special feature of the medical product plasma ONE is that it is only operated with battery power.

It does not work with or by means of mains voltage. The safe, low voltage used in battery operation provides a high level of safety, with a reduction in pathogens of up to 99.9%!

The energy charged species in cold plasma are lethal to microbes, but are completely safe for the skin. (We will be glad to provide you with our literature research for the clinical assessment of plasma ONE.

Yours sincerely


Dr./Univ. Zagreb Susanne Pavisic
Managing Director

2 Before starting use

It is essential that you read the instruction manual carefully before starting to use your new plasma ONE. As the operator of plasma ONE, strictly observe the information in the instruction manual!

We explicitly refer to the products listed by the RKI for cleaning and sterilisation and to the Medical Devices Operator Ordinance (MPBetreibV) for the use of medical products.

The plasma ONE REF 1001001/002 is suitable for use in all facilities, including domestic environments and those which are directly connected to the PUBLIC SUPPLY GRID which also supplies buildings which are used for dwelling purposes.

DIN EN 60601-1-2:2007: Section 6: EMC requirements
Emission: Limits for equipment of Class B, Group 2 according to CISPR 11
Immunity: Severity levels for none life-saving systems to be used at unshielded locations

The ambient and operating conditions are:

Ambient temperature 15 – 45° C

Relative humidity 40 – 85 %

Air pressure 800 – 1014 hPa

Operating mode: NOT FOR CONTINUOUS OPERATION

Max. activation time: 2.5 minutes

Min. deactivation time: 10 seconds

2.1 Safety information

- Note that glass is sensitive to breakage. The instrument probes must be handled with care. (Page 7-8)
- The warning information and the hygiene regulations must be communicated to the practice team by keeping records and must be complied with. (Page 9-11)
- The list of contraindications must be read carefully and must be complied with by the practice team. (Page 10)
- Use our list of recommended indications, so that you can achieve safe and successful treatment from the very start. The list of recommended indications is an accessory and can be requested at any time, free of charge.

- In case of faults or problems you can find initial information in the Solutions table.
- The service centre can be contacted under the hotline 0049(0)260393128-0 or Mail: support@plasmamedicalsystems.com.

Before using the medical product for the first time it must be disinfected or sterilised (Page 11). The regulations of the MPBetreibV must be observed.

2.2 Completeness

Please check that the contents of the package are complete.

Dental Pack

- Energy control unit
- Mains unit
- Dental converter
- Instrument probes Dental + 5a suitable sterilization basket
- Patient probe
- Recommended indications
- Instruction Manual
- Declaration of Conformity

Medical Pack

- Energy control unit
- Mains unit
- Medical converter
- Instrument probes Medical + 6a suitable sterilization basket
- Patient probe
- Recommended indications
- Instruction manual
- Declaration of Conformity

Please read the following instructions carefully before using the device.

You cannot perform treatment if the energy control unit is still connected to the mains. The treatment unit only operates in battery mode.

3 Commissioning the plasma ONE treatment unit


- 3.1 Commissioning and explanation of the energy control unit
- 3.2 Commissioning and explanation of the converter
- 3.3 Commissioning and explanation of the instrument probes
- 3.4 Commissioning and explanation of the patient probe

3.1 Commissioning and explanation of the energy control unit

Charge the energy control unit from the plug socket before the first use of the energy control unit or the complete treatment unit. To do this, plug the adapter plug connector into the

energy control unit and then plug the mains plug into the socket.

The white LED display BATTERY – CHARGING illuminates when the charging process is active! The white LED display goes out when the battery is fully charged! The green LED display POWER flashes if the charge capacity of the battery is below 30%! A fully charged battery enables at least 8 hours of active treatment.

Explanation of the display

The energy control unit is switched on with the **ON** pushbutton. The energy control unit is switched off with the **OFF** pushbutton. Press the **ON** BUTTON until the **POWER** LED display lights up green, or in case of **OFF**, until the **POWER** LED display no longer illuminates.

The levels 1 to 5 of the LED display indicate the intensity of the plasma concentration. The individual levels can be controlled with the **ON** and **OFF** switches. Click **ON** for a higher setting and **OFF** for a lower setting.

- Level 1
- Level 2
- Level 3
- Level 4
- Level 5

The recommended power levels for the particular indication is given in the list of recommended indications.

A list of indications can be found as an accessory in the packaging.

Indication	Typical instrument probe	Power Group	Time (minutes)	Application	Intervals
Gingivitis	PI04	CLASS 1	0.5 - 3	Per tooth pocket	1 - 2s, as required
	PI05				
	PI02				
Parodontitis	PI04	CLASS 1	1	Per tooth pocket	1 - 2s, as required
	PI05				
	PI02				
Pericoronitis	PI04	CLASS 1	0.5	Per tooth pocket	1 - 2s, as required
	PI05				
	PI02				
Periapicalitis	PI04	CLASS 1	1 - 1.5	Per implant	1 - 2s, as required
	PI05				
	PI02				
Periosteitis	PI04	CLASS 1	1 - 1.5	Per implant	1 - 2s, as required
	PI05				
	PI02				
Prosthetics and restorative dentistry	PI04	CLASS 1	0.5	As required	As required

The LED display **ACTIVITY** is a visual and acoustic information for the person performing treatment or the user, which indicates whether the correct and effective distance between the surface to be treated on the patient and the instrument probe has been maintained.

If the distance is correct, this LED display illuminates at intervals. If the distance is not correct, the LED display remains dark.

After 7 seconds, the treatment unit switches off for safety reasons.


The steps of the LED display **TIME/MIN** give a visual and acoustic indication of the treatment time which has elapsed. For recommended treatment times, please refer to the indications table. After 2.5 minutes, the treatment unit switches off completely for safety reasons.

The converter is started with the **START/STOP** button at the lower right of the energy control unit.

There is also a start/stop switch directly on the converter.


3.2 Commissioning and explanation of the converter

The **converter** plug connector has a push-pull lock. Plug the converter plug connector into the yellow socket at the rear of the energy control unit. The white point of the plug connector must face upwards.

3.3 Commissioning and explanation of the instrument probes

Plug the **instrument probe** (depending on the type for the corresponding indication) with the metal section into the converter. Pull the silicon ring towards the converter to secure the connection of the instrument probe. Then press **START** on the converter or on the energy control unit.

Please note that the instrument probe must be firmly inserted into the converter before pressing the **START/STOP** button on the converter or on the energy control unit!

The instrument probe type depends on the area to be treated. The wide range of plasma MEDICAL SYSTEMS® instrument probes provides access to almost any treatment point on or in the human body.

3.4 Commissioning and explanation of the patient probe

The patient probe plug connector has a push-pull lock. Plug the patient probe plug connector into the green socket at the rear of the energy control unit.

Then the patient probe is given to the patient, who holds it in his/her hand. During treatment, the patient probe creates a contact to the energy control unit, in order to monitor the indication settings.

If the patient lets go for the patient probe or does not have it in his/her hand during treatment, the treatment unit switches off after 7 seconds. This measure is for the safety of the patient and the user.

The plasma ONE is now ready for use.

Warning: Use of this cable or other accessories may impair the EMC behaviour.

In particular, the precautionary measures concerning EMC and installation must be observed (For explanatory information about the table, please refer to pages 19-21)

Explanation

The energy control unit supplies the converter with direct current. This direct current is converted into high frequencies in the converter. The high frequency which is generated is passed to the tip through the glass instrument probe. In the gap between the patient's skin and the instrument probe an electric field is created, which ionises the atoms and molecules and generates the cold plasma.

At the tip of the instrument probe an alternating electric field is set up in order to create a corona discharge (micro-discharge). The gases in the earth's atmosphere are now converted into low temperature plasmas or cold plasma.

The main gases in the earth's atmosphere are

- 78.084 % nitrogen (N₂)
- 20.946 % oxygen (O₂)
- 0.934 % argon (Ar)
- 0.036 % aerosols and trace gases

The plasma MEDICAL SYSTEMS® technology converts the earth's atmosphere, which is also referred to as the homosphere, into a disinfection medium.

4 Warning information

Only ORIGINAL plasma MEDICAL SYSTEMS® medical device components must be used.

A mains unit other than the original plasma MEDICAL SYSTEMS® mains unit is not allowed to be used.

plasma ONE may only be opened by an authorised dealer or by the manufacturer.

The plasma ONE treatment unit must only be used on patients by medical personnel.

The medical device plasma ONE is only intended for proper medical use according to the intended purpose. Do not use plasma ONE for other purposes.

Unless there are counter-indications, the discharge currents are harmless to the body. Please note the counter-indications and integrate these into your treatment concept. Discontinue treatment immediately if the patient appears to be uncomfortable.

Patients who are hypersensitive to electrical pulses may experience a tingling feeling in the treated area. However, according to present research, this is not a significant adverse side effect.

Do not switch the converter on before you have plugged the instrument probe into the converter.

Instrument probes must be cleaned and disinfected before and after use! RKI approved and recognised disinfectants and methods must be used.

Disinfection of the glass instrument probe, except for the PSKP syringe, can be carried out by autoclave sterilisation (autoclave setting 134° C, alternatively with 121° C or wiping with disinfectant) [Page 11](#).

The function of each instrument probe should be checked after cleaning and disinfection.

Handle the instrument probe with care in order to prevent damage. Glass is a sensitive material.

For example, when performing dental treatment, take care that the probe does not get caught in a periodontal pocket. For narrow or curved canals, such as root canals, the PSKP syringe should be used with Capillary Tips or Navi Tips.

Avoid any levering action with the instrument probe on the patient.

Avoid direct contact of the instrument probe with the open eye.

The unit switches off immediately in case of damage to the instrument probe. Dispose of all broken probe parts in the proper manner.

With dental treatment, sensitivity may be experienced by the patient in the case of vital teeth.

Avoid direct contact between the active instrument probe and **vital teeth** if the treatment area has **not been anaesthetised**.

First apply the probe to the area to be treated and then use the Start/Stop button to switch it on.

Do not use the treatment unit under the following conditions or in the following locations:

- If the energy control unit and/or the converter comes into contact with water or liquids.
- Where the ambient temperature and humidity exceed the values described below.

Operating conditions:

Ambient temperature

15 – 45° C

Relative humidity

40 – 85 %

Air pressure

800 – 1014 hPa.

- Where the device is subjected to vibrations or shocks.

- Do not allow the intake of liquids with the PSKP syringe during instillation. This prevents contamination of the inner core of the PSKP syringe.

- Do not use force or great pressure when treating the patient with the glass instrument probe. This prevents breakage of the glass instrument probe.

- Only use original parts from plasma MEDICAL SYSTEMS® GmbH. This prevents injury to the patient and damage to the treatment unit.

The warranty only applies for original parts. Do not use the treatment unit if electrical components are damaged and report the condition to plasma MEDICAL SYSTEMS® GmbH or the vendor or authorised dealer.

Avoid direct contact of the probe with latex materials when the device is in operation.

Some materials do not react well with low temperature plasmas. (Please note the list of plasma compatibility on [page 12](#))

Although the treatment time is relatively short, use aspiration during dental treatment, so that the unused cold plasma is eliminated from the oral cavity. This is mainly important if the duration of treatment is longer than 2 minutes.

Keep the original packaging in case you have to return the treatment unit!

4.1 Counter-indications

For the sake of procedure, possible contraindications are listed under **2, 3, 4** and **5**, although no treatment information is available, which could result in damage to health or a negative effect on the patient due to treatment with cold plasma.

Treatment must not be carried out if the patient has a pacemaker (**1**) or any other active implant.

Therefore, we recommend that patients with active implants are not given treatment!

- 1** With patients who have active implants
- 2** With patients who have severe neurological functional disorders
- 3** With patients who have severe mental functional disorders
- 4** With patients who have hypersensitivity reactions to electric pulses
- 5** With patients with severe asthma
- 6** With pregnant patients
- 7** With patients with a diagnosis of hyperthyroidism.

Please note that smaller doses should be used for the treatment of children. This primarily applies to the treatment of children's mucous membranes.

The intensity should not exceed level 2.

For all of the listed contraindications, the doctor and the patient must jointly agree whether or not plasma treatment is to be carried out.

4.2 Cleaning, disinfection & sterilisation

Cleaning and disinfection of all original parts of the plasma ONE must be carried out with disinfectants and methods approved by the RKI. During all cleaning processes, careful handling of the glass probe is essential in order to prevent breakage or damage. The metal part of the instrument probe may not be dumped down into liquid! The same applies for the electronic parts of the device and the opening in the handpiece. The whit silicon ring on the probe is used to protect the opening in the converter.

ROBERT KOCH INSTITUT

- 1** Before the intended use, all components of the plasma ONE medical device must be carefully and thoroughly cleaned with an approved and entirely manual method.
- 2** After this, the components of the device are disinfected with an RKI-approved agent for wiping disinfection (e. g. isopropanol or ethanol) or an approved thermal method at 121° C.
- 3** Sterilisation of the glass probes should only be carried out for use of the plasma ONE in a critical area of use (e.g. surgery or dentistry), using the sterilisation method approved and required by the RKI. Sterilisation of all glass instrument probes, except for the PSKP, should be performed in an autoclave at 134° C (up to 100x sterilisation cycles).

After each sterilisation and disinfection, the instrument probes must be dried and properly stored. Their function must be visually inspected before the next use.

All parts of the device completely fulfil the safety requirements. Because of this, the device must not be used with instrument probes or other accessories which are not included in the scope of delivery or which have not been tested by the manufacturer.

Warning:
According to DIN IEC 60601-1 the plasma ONE treatment unit must not be opened or structurally modified without the consent of the manufacturer!

4.3 Service & maintenance

The device is maintenance free. However, the European Ordinance and the Medical Products Act MPG for medical devices require a safety inspection every two years.

The plasma MEDICAL SYSTEMS® products must be inspected in an authorised service center or by the manufacturer according to the MPBetreibV. In order to prevent possible damage, please do not use any components or device parts which are not plasma MEDICAL SYSTEMS® original parts.

Please use the original packaging or the case for storage or transport of the plasma MEDICAL SYSTEMS® products.

The performance of glass instrument probes will decrease in the course of time, especially if the instrument probes are frequently sterilised. Because of this, we recommend that the instrument probes are replaced every 12 months, depending on how often they are used and sterilised.

Warning information

Information for users

The MEDICAL SYSTEMS® GmbH warranty does not apply to faults and breakages which are due to improper cleaning, disinfection or sterilisation.

The white silicon ring on the instrument probes (except for the PSKP syringe) is used to protect the opening on the converter. It can be removed for sterilisation purposes.

Disposal of the electronic components

(Applicable in the European Union and other European states with waste sorting systems.)

This symbol on a product or its packaging means that the product may not be disposed of as domestic waste. Rather, it should be taken to a collection point for the recycling of electronic scrap. By disposing of this product correctly, you can help to prevent environmental damage and damage to human health which could

occur due to the improper disposal of this product. Recycling of materials also helps to conserve natural resources. The energy control unit must not be opened if the device is defective.

The service life of the plasma ONE is declared to be at least 6 years, providing that it is properly handled and used according to the rules which are described in the instruction manual.

Warning: The device must not be used in the immediate vicinity of other devices. Portable and mobile HF communication equipment may influence medical electrical devices.

4.4 List of plasma compatibility

Assessment of chemical effects

- 1 **Very good** = No effect
- 2 **Good** = Slight effect
Slight corrosion or discolouration
- 3 **Medium** = Moderate effect
Not recommended for continuous use. Softening, loss of strength or swelling may occur.
- 5 **Poor** = Not recommended

Assessment of materials

- 1 **Very good**
Butyl, Chemraz, CPVC, Durachlor-51, Durlon 9000, EPDM, EPR, ethylen propylene, fluorsilicone, galvanised steel, water in the air, glass, Hastelloy-C®, Hypalon®, Hytrel®, Inconel, Kalrez, Kel-F® (PCTFE, Norprene, PEEK, polycarbonates, polyurethane, rollable, PTFE (Teflon®, PVDF (Kynar®), Santoprene, silicons, stainless steel – 316, titanium, Tygon®, Vamac, Viton®
- 2 **Good**
Aluminium, copper, bronze, LPDE, polysulfide, PVC, polyacrylate, stainless steel – 304
- 3 **Medium**
Acetal (Delrin®), cast iron, galvanised steel, water in glass, monel metal, neoprene, polypropylene
- 5 **Poor**
Buna-N (Nitrile) > Severe effect, magnesium, natural rubber > Severe effect, Nylon > Severe effect, steel (construction steel, HSLA), zinc

5 Information for users

5.1 Successful treatment – step-by-step

- 1 Place the converter and the patient probe in the relevant holder on the instrument table.
- 2 Check that the battery is charged and that the mains plug is not in the socket.
- 3 Decide which instrument probe you wish to use for the treatment. Plug the instrument probe with the metal piece onto the converter and pull the silicon ring up to the edge of the converter.
- 4 Place the patient probe in the patient's hand.

With the patient probe, the patient is actively involved in the success of the treatment.

Explain to the patient that slight pulses of current which have a control function will be conducted through his/her hand to the energy control unit.

Because each patient has an individual sensitivity a measurement is made, which monitors the pre-set power criteria.
If there is any change during treatment, the **ACTIVITY** LED display switches off immediately.

The person performing the treatment can correct this immediately and change the distance to the treatment area or reduce the power level.

Changes may occur if the distance between the instrument probe and the treated area is not correct, or the patient suffers from hypersensitivity or sensitivity.

If the patient removes the patient probe from his/her hand, he/she has a direct effect on the discontinuation of treatment.

- 5 Switch on the treatment unit at the converter or the energy control unit.
- 6 Set the power level from 1 to 5 with the **ON** button. Click the **ON** button until you have reached the required level.
- 7 Then apply the tip of the instrument probe to the treatment area.
Take care to maintain a distance of up to 5 mm between the treatment area and the tip of the probe.
Examples of how the various instrument probes can be used on the treatment area in dermatology (including for tattoos).
- 8 When you have completed the treatment, proceed as follows:
Switch off the converter at the converter or press the **START/STOP** button on the lower right of the energy control unit.

- 9 Carefully pull the instrument probe out of the converter. Put down the instrument probe in a safe place. **Disinfect the used instrument probe before using it on a new patient.**

Information for users

Information for users

If you continue treatment with a different type of instrument probe, proceed as described above.

Plug the instrument probe with the metal part into the converter with a slight rotation and start treatment by pressing the **START/STOP** button on the energy control unit or directly on the converter.

Instrument probes may only be exchanged when the converter is switched off. Instrument probes must be sterilised before use

The power levels and treatment times, which are recommended in the indications table, result from an extensive in-vitro study with the relevant pathogens, bacteria and fungi, as well as from clinical experience from daily work in practice.

The success of the treatment depends on the consistent guidance of the instrument probe. The gap between the treatment area and the tip of the instrument probe should be about less than 2 mm.

The condition of the patient described above, or his/her individual sensitivity and physical reaction can change the pre-set intensity. These changes are recorded by the use of the patient probe. They are communicated to the energy control unit and displayed as described.

The person performing the treatment recognises that he/she needs to improve the placement of the instrument probe.

It is beneficial if the person performing the treatment or the user also creates a change of air by means of a slight movement. A rigid position is just as unfavourable as too large a gap between the treated area and the instrument probe.

If the technological process detects an optimal handling, the gas in the glass instrument probe illuminates in a high-tech orange and transparently indicates the correct functioning to all concerned.

5.2 Experience in practice

The individual condition of the patient determines the intensity!

In case of sensitive patients or parts of the body, a combination of levels 1 to 2 and a maximum treatment time of 2 minutes is sufficient.

We recommend this setting for children.

During treatment, the patient may experience a slight tingling sensation. This feels as though a wound is drying in air and causes slight itching. As with natural wound healing, this is a good and positive reaction.

According to current research there is no danger of an overdose of cold plasma, as low temperature plasmas decompose within 20 seconds.

Overview of intensities or power levels

The intensity settings or power levels are divided into three main groups:

- **Group III High intensity**
Level 4 and 5
- **Group II Medium intensity**
Level 3 to 4
- **Group I Low intensity**
Level 1 to 2

Group II Medium intensity

Level 3, up to or/and level 4. Used for improved **healing and regeneration**.

If the patient is particularly sensitive, then Group III is indicated at level 3 with a longer treatment time!

- Aphthae
- Burns
- Wound healing
- After surgical operations
- Pressure sores
- Haematoma

Group I Low intensity

Level 1 and/or level 2

Used for sensitive patients and hypersensitive treatment areas.

- Hypersensitive teeth
- Patients with sensitive tissue
- Prophylaxis
- Children

Detailed explanations

Group III High intensity

Level 4 and level 5 illuminate together. This is used for example if a high level of **disinfection** is required and for a high tissue stimulation.

- Abscesses
- Bacterial infections
- Ulcers
- Virus infections, e.g. Herpes simplex
- Candida albicans (fungal infections)
- Wound infection, especially poorly healing wounds after secondary infection
- Long-term bleeding
- Fistulas
- Gangrene
- Gingivitis
- Granuloma
- Implantation
- Periodontitis
- Periimplantitis
- Disinfection of the pulp
- Root canal treatment
- Post and crown disinfection
- Caries

As a matter of principle, it should be observed whether the patient indicates an unpleasant feeling! Hypersensitivity should be compensated for with a lower level/intensity and a longer treatment time.

Patient reports

Inform each patient of the result of the treatment by means of a treatment report. This makes the versatile effect of the treatment transparent. This helps both, you and your patients to establish the convenient, pain-free treatment with cold plasma as a routine standard treatment. This results in satisfied patients and an accepted criterion for invoicing private medical treatment.

A brief overview of the various instrument probes:

Dermatology, gynaecology, surgery

A recommendation of which type of instrument probe is suitable for which area of use can be found in the accessory list of recommended indications.

Dentistry


Installation of the PSKP syringe

1 Before use, fill the PSKP syringe (E4) with air. To do this, hold the PSKP syringe firmly with one hand and with the other hand pull the lower section with the metal end (E4a) downwards. The air inside the transparent air chamber is later converted into cold plasma.

2 The NaviTips (E4c) or Capillary (E4b) Tips must be pushed onto the PSKP syringe (E4) before the PSKP probe is plugged into the converter.

Plug the PSKP syringe with the metal end (E4a) into the converter (D3).

- Now switch on the energy control unit by pressing the **ON** button until the POWER LED display illuminates green.
- Via the ON button you can reach level 5.
- Then switch on the converter with the **START/STOP** button.
- Hold the top section of the PSKP syringe (E4) with two fingers for approx. 10 seconds. This action accelerates the cold plasma conversion process in the air chamber. The patient probe has to be maintained during these 10 seconds in user's hand.
- After approx. 10 seconds, switch off the converter at the converter or the energy control unit, because the filling process is now complete.
- **Now remove the PSKP syringe from the converter and inject it into the relevant treatment area.**

Using the PSKP syringe

The versatility of the various instrument probes will only be used in dental medicine.

6.1 Solution table

Problem	Überprüfung	Lösung
The instrument probe is not active, although the treatment unit is switched on.	<p>Was the Start/Stop button pressed for at least 3 seconds?</p> <p>Is the converter connected to the energy control unit?</p> <p>Is the instrument probe damaged in any way? If yes...</p> <p>Is the patient holding the patient probe during treatment?</p>	<p>Press the START/STOP button on the energy control unit or directly on the converter until the power indicator light illuminates!</p> <p>Check whether the converter plug at the rear of the energy control unit is properly plugged into the socket! The arrow on the plug must be visible at the top!</p> <p>Replace the defective instrument probe!</p> <p>Please plug the plug connector of the patient probe into the rear of the energy unit and place the probe in the patient's hand.</p>
The treatment unit switches off if the instrument probe comes into contact with the patient...	<p>Has the glass of the instrument probe been damaged during the treatment?</p> <p>Patients can store more current or can generate a greater electronic charge than has been preprogrammed as standard.</p> <p>This is reported to the energy control unit by the patient probe.</p> <p>Did the end of the instrument probe which is plugged into the converter become wet during the treatment?</p>	<p>Remove the broken glass from the patient and dispose of the defective instrument probe in a proper manner!</p> <p>Restart treatment by setting the intensity 1 to 2 levels lower and extend the treatment time!</p> <p>Dry the entire instrument probe, especially the metal piece!</p>
The treatment unit switches off automatically after charging of the PSKP...	<p>Is there liquid or other material in the probe?</p> <p>Has the probe been properly plugged into the handpiece?</p>	<p>Replace the PSKP with a new syringe!</p> <p>Reinstall the PSKP!</p>

7.1 Explanation of symbols

Observe the instruction manual

Mains unit, protection class II
When in operation, the treatment unit only has an internal power supply with a safe low voltage.

Use without earthing,
application unit type BF

C E Declaration of Conformity

These devices comply with the European Directive 2002/96/EU for (old) electrical devices (disposal of electrical and electronic scrap)

For identification of the positive and negative contacts (polarity) of a DC supply or the positive and negative contacts of a device which can be supplied with a DC voltage.

For indication of an electrode holder with a button, e.g. for storage, use and connection.

IP 20

IP Code 20
1. Figure (2)
Protected against solid foreign bodies with a diameter larger than 12.5 mm
2. Figure (0)
No protection against water

The HF transmitters contain or use electromagnetic HF energy for diagnosis or treatment: Labeling with the symbol according to EN 60417-5140 for non-ionising radiation

7.2 Operating mode

Operating mode plasma ONE medical device

Li-ion	battery pack
	7.2 V/2.6 Ah
Power output	Maximum 5 W
Device fuse	1 A Polymer, built-in/ self-resetting
	30R110 / E183209
Repetition frequency	420 Hz – 1220 Hz
Converter pulse width	2.5 µs – 10 µs
Maximum patient current DENTAL	< 100 µA
Maximum patient current MEDICAL	< 100 µA
Automatic power switch-off	t < 0.1 s. max 200 µA
Maximum operating altitude	< 2000 m
Protection class according to EN 60529	IP 20
606201-1 8.9.1.9	300 A – 1500

Operating mode Charging with mains unit on the medical device

Mains unit input	100 – 240 V AC
Insulation voltage	corresponds to standard EN60950/EN60601
Mains unit output	12 V DC

Dimensions and weight

Length x width x max. height	22 cm x 19 cm x 8 cm
Weight	1,250 gr

Ambient and storage conditions

Ambient temperature	15 – 45° C
Relative humidity	40 – 85 %
Air pressure	800 – 1014 hPa

Storage and transport conditions

Ambient temperature	5 – 60° C
Relative humidity	20 – 85 %
Air pressure	800 – 1014 hPa
Dispatch/Transport	only in original case from plasma MEDICAL SYSTEMS® GmbH

7.3 EMC information

Directives and manufacturer's declaration – Electromagnetic interference		
The plasma ONE REF 1001001/002 is intended for operation in one of the ELECTROMAGNETIC ENVIRONMENTS listed below. The customer or the user of the plasma ONE REF 1001001/002 must ensure that it is operated in such an environment.		
Radio interference measurements	Compliance	Electromagnetic environment directives
HF emissions as per CISPR 11	Group 2	The plasma ONE REF 1001001/002 must emit electromagnetic energy in order to fulfil its intended function. This may affect adjacent electronic devices.
HF emissions as per CISPR 11	Class B	The plasma ONE REF 1001001/002 is suitable for use in all facilities, including domestic environments and those which are directly connected to the PUBLIC SUPPLY GRID which also supplies buildings which are used for dwelling purposes.
Emission of harmonics as per IEC 61000-3-2	Class A	
Emission of voltage fluctuations /flicker as per IEC 61000-3-3	Compliant	

Recommended safety distances between portable and mobile HF telecommunication devices and the plasma ONE REF 1001001/002			
The plasma ONE REF 1001001/002 is intended for use in an electromagnetic environment in which HF interference is under control. The customer or user of the plasma ONE REF 1001001/002 can help to reduce electromagnetic interference by complying with the minimum distance between portable and mobile telecommunication devices (transmitters) and the plasma ONE REF 1001001/002 – depending on the output power of the communication device as stated below.			
Rated power of transmitter W	Safe distance depending on transmission frequency m		
	150 kHz to 80 MHz d = 1.2√P	80 MHz to 800 MHz d = 1.2√P	800 MHz to 2.5 GHz d = 2.3√P
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters whose rated power is not stated in the table above, the recommended safe distance d in metres (m) can be calculated with the equation for the relevant column, whereby P is the maximum rated power of the transmitter in Watt (W) according to the information from the manufacturer. NOTE 1 With 80 MHz and 800 Hz, the higher frequency range applies NOTE 2: These guidelines may not be applicable in all cases. The dissipation of electromagnetic effects is influenced by absorption and reflection by the building, objects and people.			

ONE REF 1001001 = Medical
ONE REF 1001002 = Dental

7.3 EMC information

Directives and manufacturer's declaration - Resistance to electromagnetic interference			
The plasma ONE REF 1001001/002 is intended for operation in one of the ELECTROMAGNETIC ENVIRONMENTS listed below. The customer or the user of the plasma ONE REF 1001001/002 must ensure that it is operated in such an environment.			
Interference immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment directives
DISCHARGE OF STATIC ELECTRICITY (ESD) as per IEC 61000-4-2	±6 kV Contact discharge ±8 kV Air discharge	±6 kV Contact discharge ±8 kV Air discharge	Floors should consist of wood or concrete or be covered with ceramic tiles. If the floor is covered with synthetic materials, the relative humidity be at least 30%.
Rapid transient electric interference/Bursts as per IEC 61000-4-4	± 2 kV for mains cables ± 1 kV for input and output cables	± 2 kV for mains cables Not applicable	The quality of the supply voltage should correspond to that of a typical commercial or hospital environment.
Peak voltages/Surges as per IEC 61000-4-5	± 1 kV voltage, external conductor - external conductor ±2 kV External conductor-earth	± 1 kV voltage, external conductor - external conductor Not applicable	The quality of the supply voltage should correspond to that of a typical commercial or hospital environment.
Voltage drops, brief interruptions and fluctuations of the supply voltage as per IEC 61000-4-11	< 5% U_T (> 95% breakdown of U_T) for ½ period 40% U_T (60% breakdown of U_T) for 5 periods 70% U_T (30% breakdown of U_T) for 25 periods < 5% U_T (> 95% breakdown of U_T) for 5 seconds	Not applicable	Not applicable
Magnetic fields with the supply frequency (50/60 Hz) as per IEC 61000-4-8	3 A/m	3 A/m	Magnetic fields with the supply frequency should correspond to the typical values for a commercial and hospital environment.
NOTE: U_T is the alternating mains voltage before application of the test level.			

ONE REF 1001001 = Medical
ONE REF 1001002 = Dental

7.3 EMC information

Directives and manufacturer's declaration - Resistance to electromagnetic interference			
The plasma ONE REF 1001001/002 is intended for operation in one of the ELECTROMAGNETIC ENVIRONMENTS listed below. The customer or the user of the plasma ONE REF 1001001/002 must ensure that it is operated in such an environment.			
Interference immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment directives
			Portable and mobile radio devices should not be used at a closer distance to the plasma ONE REF 1001001/002, including the cables, than the recommended protection distance which is calculated according to the relevant equation for the transmission frequency. Recommended safe distance:
Conducted HF interference as per IEC 61000-4-6	3 V_{eff} 150 kHz to 80 MHz	3 V_{eff}	$d = 1.2\sqrt{P}$
Radiated HF interference as per IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2.5 GHz
			Where P is the rated power of the transmitter in Watt (W) according to the details from the manufacturer, and d is the recommended safe distance in metres (m). For all frequencies, the field strength of fixed radio transmitters should be less than the compliance level ^b as level as determined by local investigation ^a . Interference is possible in the vicinity of devices labelled with the following symbol. 
NOTE 1 With 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2: These guidelines may not be applicable in all cases. The dissipation of electromagnetic effects is influenced by absorption and reflection by the building, objects and people.			
a The field strength of fixed transmitters, e.g. the base stations of radio telephones, and mobile terrestrial radio devices, amateur radio stations, AM and FM radio and television transmitters cannot be theoretically determined in advance. In order to determine the electromagnetic environment with regard to fixed radio transmitters, a study of the location should be considered. If the measured field strength at the location in which the plasma ONE REF 1001001/002 is to be used exceeds the compliance levels stated above, the device should be observed in order to provide proof of correct functioning. If unusual performance characteristics are observed, additional measures may be necessary, e.g. a change of orientation or a different location of the plasma ONE REF 1001001/002. b Over the frequency range of 150 kHz to 80 MHz the field strength should be less than 3 V/m.			

8.1 ISO 13485, ISO 9001 and EC-Certificate

Warranty

plasma MEDICAL SYSTEMS® GmbH
 plasma ONE treatment unit

Energy control unit	24 months
Converter	24 months
Instrument probe	12 months
Patient probe	12 months

Faults and breakages due to misuse are not covered by this warranty.

Faults and breakages which are due to improper disinfection or sterilisation processes are also not covered by this warranty.

The warranty is only valid if the plasma MEDICAL SYSTEMS® treatment unit is serviced by the dealer's authorised personnel or by the manufacturer.

The warranty is void if accessory components other than those which are originally supplied or tested by the manufacturer are used, or if the treatment unit is opened without authorisation.

The instructions in this instruction manual must be observed. Otherwise the warranty becomes void.

plasma MEDICAL SYSTEMS® GmbH

Viktoriaallee 1 | D-56130 Bad Ems/Germany

Phone +49 (0) 2603-93128-0

Fax +49 (0) 2603-93128-28

info@plasmamedicalsystems.com

www.plasmamedicalsystems.com

01/2013-07/DE/EN

© Copyright plasma MEDICAL SYSTEMS® GmbH

CE
0535

German Medical Devices Act (MPG) in the version
of 19 October 2012 (Federal Law Gazette I p. 2192)

The EC directive 93/42/EEC
in particular Annex I and Annex V.
Applied harmonised standards
DIN EN ISO 13485:2012
DIN EN ISO 9001:2008
IEC 60601-1, 3rd Edition