







Preis trifft Leistung.

PERFECT MATCH!

Die neue Endfeilen-Generation von EdgeEndo









EDGEONE
FIRE™

Technische Übersicht

Reziprok 30° Rechtslauf und 150° Linkslauf	Spitzengröße	Max. Durchmesser	16mm	13mm	9mm	6mm	3mm	1mm	Drehmoment g/cm Ncm
 Klein	0.20	1.20	1.20	0.87	0.72	0.54	0.37	0.25	408 – 530 g/cm 4,0 – 5,2 N/cm
 Normal	0.25	1.20	1.20	0.98	0.82	0.61	0.43	0.30	408 – 530 g/cm 4,0 – 5,2 N/cm
 Mittel	0.35	1.20	1.20	0.95	0.82	0.63	0.49	0.40	408 – 530 g/cm 4,0 – 5,2 N/cm
 Groß	0.45	1.20	1.20	0.92	0.82	0.66	0.56	0.50	408 – 530 g/cm 4,0 – 5,2 N/cm

EDGETAPER **EDGETAPER**
PLATINUM™

Technische Übersicht

	Spitzengröße	Max. Durchmesser	16mm	13mm	9mm	6mm	3mm	1mm	Drehmoment (Ncm)	Drehzahl (U/Min)
 SX	0.19	1.20	-	1.20	1.15	0.78	0.40	0.24	2,5 - 3,0	300 - 400
 S1	0.17	1.20	1.20	1.02	0.68	0.49	0.30	0.21	2,5 - 3,0	300 - 400
 S2	0.20	1.20	1.20	1.01	0.77	0.57	0.37	0.25	2,5 - 3,0	300 - 400
 F1	0.20	1.20	1.14	0.98	0.77	0.59	0.42	0.28	2,5 - 3,0	300 - 400
 F2	0.25	1.20	1.15	1.01	0.81	0.65	0.49	0.35	2,5 - 3,0	300 - 400
 F3	0.30	1.20	1.18	1.06	0.89	0.73	0.57	0.41	2,5 - 3,0	300 - 400
 F4	0.40	1.20	1.20	1.06	0.87	0.72	0.57	0.46	2,5 - 3,0	300 - 400
 F5	0.50	1.20	1.20	1.06	0.87	0.75	0.63	0.55	2,5 - 3,0	300 - 400










Preis trifft Leistung.

PERFECT MATCH!

Die neue Endofeilen-Generation von EdgeEndo








EDGEFILE® X7.04

Technische Übersicht

	Spitzengröße	Taper	Max. Durchmesser	16mm	13mm	9mm	6mm	3mm	1mm	Drehmoment Ncm	Drehzahl (U/Min)
	17	0.04	0.8	0.8	0.69	0.53	0.41	0.29	0.21	2,45-4,02	300-500
	20	0.04	0.8	0.8	0.72	0.56	0.44	0.32	0.24	2,45-4,02	300-500
	25	0.04	0.8	0.8	0.77	0.61	0.49	0.37	0.29	2,45-4,02	300-500
	30	0.04	0.8	0.8	0.8	0.66	0.54	0.42	0.34	2,45-4,02	300-500
	35	0.04	0.8	0.8	0.8	0.71	0.59	0.47	0.39	2,45-4,02	300-500
	40	0.04	1.0	1.0	0.92	0.76	0.64	0.52	0.44	2,45-4,02	300-500
	45	0.04	1.0	1.0	0.97	0.81	0.69	0.57	0.49	2,45-4,02	300-500

EDGEFILE® X7.06

Technische Übersicht

	Spitzengröße	Taper	Max. Durchmesser	16mm	13mm	9mm	6mm	3mm	1mm	Drehmoment Ncm	Drehzahl (U/Min)
	17	0.06	1.0	1.0	0.95	0.71	0.53	0.35	0.23	2,45-4,02	300-500
	20	0.06	1.0	1.0	0.98	0.74	0.56	0.38	0.26	2,45-4,02	300-500
	25	0.06	1.0	1.0	1.0	0.79	0.61	0.43	0.31	2,45-4,02	300-500
	30	0.06	1.0	1.0	1.0	0.84	0.66	0.48	0.36	2,45-4,02	300-500
	35	0.06	1.0	1.0	1.0	0.89	0.71	0.53	0.41	2,45-4,02	300-500
	40	0.06	1.0	1.0	1.0	0.94	0.76	0.58	0.46	2,45-4,02	300-500
	45	0.06	1.0	1.0	1.0	0.99	0.81	0.63	0.51	2,45-4,02	300-500